



## ADVANCED FILL OVERVIEW

	◀	<b>HYPOTHERMIA BLANKET</b>		▶	<b>COMPRESSED GAUZE</b>		▶	<b>PRESSURE DRESSING</b>
	▶	<b>FLAT FOLD TAPE</b>		▶	<b>TRAUMA GLOVES</b>		▶	<b>DECOMPRESSION NEEDLE</b>
	▶	<b>CHEST SEAL</b>		▶	<b>QUIKLOT HEMOSTATIC DRESSING</b>		▶	<b>NASOPHARYNGEAL AIRWAY</b>

**When folding excess medical component packaging, take care to not fold, crease or crush the product inside. Please contact Customer Service if you have questions.**



**1) Insert Placement**  
To begin, place the inner Trauma Kit out with the pull tab (medical insignia) facing as shown. Then insert Medical items A-C in alphabetical order.



**2) Medical Items Placement**  
Stack items as shown below

- A** CHEST SEAL
- PRESSURE DRESSING
- QUIKLOT HEMOSTATIC DRESSING
- B** FLAT FOLD TAPE
- HYPOTHERMIA BLANKET
- COMPRESSED GAUZE
- C** DECOMPRESSION NEEDLE
- D** NASOPHARYNGEAL AIRWAY
- E** TRAUMA GLOVES

**REQUIRES FOLDING**

**Securing Your Medical Items**

1. Fold Trauma Kit inner pouch
2. Insert into Pouch
3. Secure with Hook and Loop Flap





## PRO FILL OVERVIEW



◆ HYPOTHERMIA BLANKET



◆ COMPRESSED GAUZE



◆ PRESSURE DRESSING



◆ FLAT FOLD TAPE



◆ TRAUMA GLOVES



◆ CHEST SEAL



◆ QUIKCLOT HEMOSTATIC DRESSING



**1) Insert Placement**

To begin, place the inner Trauma Kit out with the pull tab (medical insignia) facing as shown. Then insert Medical items A-C in alphabetical order.



**2) Medical Items Placement**

Stack items as shown below



**REQUIRES FOLDING**

**Securing Your Medical Items**

1. Fold Trauma Kit inner pouch
2. Insert into Pouch
3. Secure with Hook and Loop Flap



**When folding excess medical component packaging, take care to not fold, crease or crush the product inside. Please contact Customer Service if you have questions.**



## ESSENTIALS FILL OVERVIEW



▶ HYPOTHERMIA BLANKET



▶ COMPRESSED GAUZE



▶ FLAT FOLD TAPE



▶ TRAUMA GLOVES



▶ PRESSURE DRESSING



▶ QUIKCLOT HEMOSTATIC DRESSING


**1) Insert Placement**

To begin, place the inner Trauma Kit out with the pull tab (medical insignia) facing as shown. Then insert Medical items A-C in alphabetical order.



**2) Medical Items Placement**

Stack items as shown below



**REQUIRES FOLDING**

**Securing Your Medical Items**

1. Fold Trauma Kit inner pouch
2. Insert into Pouch
3. Secure with Hook and Loop Flap



**When folding excess medical component packaging, take care to not fold, crease or crush the product inside. Please contact Customer Service if you have questions.**